



## TEST REPORT NO 125321/23/GDY/2

SFD SPÓŁKA AKCYJNA GŁOGOWSKA 41 45315 OPOLE		Sample (according to declaration of Client) Sample description: ALLNUTRITION MAGNESIUM 5 FORMS + B6 (P-5-P) 100 kaps Batch: AN230101 Expiry date: 10.01.2025		
Sample reception date:	10.03.2023	Sample status: no objections		
Start of analysis 15.03.2023		1		
End of analysis	20.03.2023	Complex received from the Olient	mple received from the Client	
Test report date 21.03.2023		Sample received from the Client		
Test Method		Unit	Result	
* Number of yeasts and moulds PN-ISO 21527-2:2009 (withdr				
Number of yeasts		cfu/g	<1,0x101	
Number of moulds		cfu/g	<1,0x101	
* Presence of coagulase-positive staphylococci (Staphylococcus aureus and other species) in 1 g PN-EN ISO 6888-3:2004; PN-EN ISO 6888-3:2004/AC:2005		in 1 g	Not detected	
* Presence of Escherichia coli in 1 g PN-ISO 7251:2006		in 1 g	Not detected	
* Presence of Salmonella spp. in 25 g PN-EN ISO 6579-1:2017-04; PN-EN ISO 6579-1:2017-04/A1:2020-09		in 25 g	Not detected	
* Presence of Listeria monocytogenes in 25 g PN-EN ISO 11290-1:2017-07		in 25 g	Not detected	
* Pyrrolizidine alkaloids <sup>3) 4)</sup> PB-498 ed. I of 23.05.2022				
Echimidine		μg/kg	< 5,0 (5,0 ± 1,8)	
Echimidine N-oxide		μg/kg	< 5,0 (5,0 ± 1,8)	
Echinatine N-oxide		μg/kg	< 5,0 (5,0 ± 1,8)	
Europine		μg/kg	< 5,0 (5,0 ± 1,8)	
Europine N-oxide		μg/kg	< 5,0 (5,0 ± 1,8)	
Heliosupine		μg/kg	< 5,0 (5,0 ± 1,8)	
Heliosupine N-oxide		μg/kg	< 5,0 (5,0 ± 1,8)	
Heliotrine		μg/kg	< 5,0 (5,0 ± 1,8)	
Heliotrine N-oxide		μg/kg	< 5,0 (5,0 ± 1,8)	
Intermedine		μg/kg	< 5,0 (5,0 ± 1,8)	
Intermedine-N-oxide (sum of intermedine-N-oxide and indicine-N-oxide as intermedine-N-oxide)		μg/kg	< 5,0 (5,0 ± 1,8)	
Lasiocarpine		μg/kg	< 5,0 (5,0 ± 1,8)	
Lasiocarpine N-oxide		μg/kg	< 5,0 (5,0 ± 1,8)	





## **TEST REPORT NO 125321/23/GDY/2**

Lycopsamine (sum of lycopsamine, indicine and echinatine as lycopsamine)	μg/kg	< 5,0 (5,0 ± 1,8)
Lycopamine N-oxide	μg/kg	< 5,0 (5,0 ± 1,8)
Retrorsine (sum of retrorsine and usaramine as retrorsine)	μg/kg	< 5,0 (5,0 ± 1,8)
Retrorsine-N-oxide	μg/kg	< 5,0 (5,0 ± 1,8)
Rinderine	μg/kg	< 5,0 (5,0 ± 1,8)
Rinderine N-oxide	μg/kg	< 5,0 (5,0 ± 1,8)
Senecionine	μg/kg	< 5,0 (5,0 ± 1,8)
Senecionine-N-oxide (sum of senecionine-N-oxide and integerrimine-N -oxide as senecionine-N-oxide)	μg/kg	< 5,0 (5,0 ± 1,8)
Seneciphylline (sum of seneciphylline and spartioidine as seneciphylline)	μg/kg	< 5,0 (5,0 ± 1,8)
Seneciphylline-N-oxide (sum of seneciphylline-N-oxide and spartioidine N-oxide as seneciphylline-N-oxide)	μg/kg	< 5,0 (5,0 ± 1,8)
Senecivernine (sum of senecivernine and integerrimine as senecivernine)	μg/kg	< 5,0 (5,0 ± 1,8)
Senecivernine N-oxide	μg/kg	< 5,0 (5,0 ± 1,8)
Senkirkine	μg/kg	< 5,0 (5,0 ± 1,8)
Usaramine N-oxide	μg/kg	< 5,0 (5,0 ± 1,8)
Sum of pyrrolizidine alkaloids	μg/kg	below quantification limit
* Vitamin B6 <sup>1)</sup> PN-EN 14164:2014-08		
Vitamin B <sub>6</sub> (pyridoxine hydrochloride)	mg/capsule	1,82
Vitamin B <sub>6</sub> (pyridoxine)	mg/capsule	1,50
* Magnesium (Mg) ²) PB-36/ICP ed. VII of 07.07.2022	mg/capsule	105

1) Capsule weight declared by the Client: 900 mg.

<sup>2)</sup> The weight of capsules declared by the Client: 900 mg.

3) The lower limit of the measuring range of the accredited method, which is also the limit of quantification set by the Laboratory.

<sup>4)</sup> Limit of quantification 5,0 (5,0 ± 1,8) μg/kg.

Authorized by: Anna Polanin, Manager, Microbiology Laboratory Ewa Ostrach-Grzybowska, Analysis Expert, Vitamin Analysis Laboratory Katarzyna Szpinda, Analysis Expert, Spectrometry Laboratory Małgorzata Grzegorek, Analysis Expert, Liquid Chromatography Laboratory

The test report bears the certified electronic seal of J.S. Hamilton Poland Sp. z o.o. Laboratory address: Chwaszczyńska 180, 81-571 Gdynia Ks. Stanisława Kujota 8, 70-605 Szczecin

## THE END OF THE REPORT

The results refer only to the samples received. When a measurement uncertainty is given, it is an expanded uncertainty estimated for a coverage factor k=2 at 95% confidence level and is not including sampling uncertainty, unless otherwise stated. When the conformity is stated J.S. Hamilton Poland Sp. z o.o. applies the simple acceptance decision rule in accordance with ILAC-G8:09/2019, unless otherwise reported. If the "result" column of the accredited method contains a record: "<" or ">", it means, that it is the test outcome directly related to the lower or upper limit of the measuring range of the accredited method respectively. In such a case, the Laboratory presents the opinion and interpretation in the "statement of conformity" column, which is based on the obtained test outcome. This test report may not be copied in part without the prior written permission of J.S. Hamilton Poland Sp. z o.o. The responsibility of J.S. Hamilton Poland Sp. z o.o. does not permit the use of the PCA accredited method symbol. AB 079 by customers, subcontractors, external service providers and other third parties. For further information please refer to the PCA document - DA-02. The service confirmed by this report is subject to the General Terms and Conditions of Services of J.S. Hamilton Poland Sp. z o.o. published on www.hamilton.com.pl.

Page 2/3





## TEST REPORT NO 125321/23/GDY/2

\* Test method accredited # Test performed by external provider

Page 3/3